

INTERNAL USE ONLY
— DO NOT PRINT THIS PAGE —

MMF-RM7318 Rev 1 - Package Insert, IFU, HemCon OneStop Bandage
MS-RM4018
Size: 4.25" x 5.5", folded to 4.25" x 2.75"
Color: Black, PMS 185

HEMCON[®]

OneStop[™]
Bandage

The blood stops here.[™]

Instructions for Use:

1. Apply dressing directly over actively bleeding wound with backing away from wound. If necessary, the bandage can be cut or folded to the size of the wound.
2. Hold pressure until bleeding is controlled.
3. Apply outer bandage wrap (not included) to secure dressing on wound site, if required.
4. Remove dressing within 48 hours with saline or sterile water.

Intended Use/Indications: The OneStop™ Bandage is intended for the local management of bleeding wounds and to provide a barrier to bacterial penetration of the dressing in all patients and for the promotion of rapid control (hemostasis) of bleeding in patients following hemodialysis and for those on anticoagulation therapy. The dressing is indicated for the following wounds: lacerations, abrasions, nose-bleeds, surgical debridement sites, skin surface puncture sites, vascular procedure sites, and sites involving percutaneous catheters, tubes, and pins.

Description: The OneStop™ Bandage is a wound dressing manufactured from chitosan. When applied directly to the wound, the OneStop™ Bandage controls bleeding. The OneStop™ Bandage offers an antibacterial barrier against a wide range of gram-positive and gram-negative organisms, including antibiotic resistant *Staphylococcus aureus* (MRSA), *Enterococcus faecalis* (VRE), and *Acinetobacter baumannii*. Only single strains of most species mentioned have been studied.

Cautions/Warnings: Do not use if seal is broken and/or contents are wet. Do not apply over eyes. Use promptly after opening. Contains chitosan from shellfish. For external temporary use only, do not implant. Not for consumption, do not eat. Single use only. Discard any unused portion of the product. Do not resterilize. Not made with natural rubber latex.

Seek Medical Attention: If bleeding does not stop or redness or swelling occurs. If you suspect an allergic reaction. For severe wounds that have breached the dermis, to avoid potential complications due to internal bleeding.

Storage Conditions: Best stored at room temperature.



Manufactured for:
Tricol Biomedical, Inc.

720 SW Washington Street, Suite 200,
Portland, OR 97205 U.S.A.

+1.503.245.0459 | Toll Free 877.247.0196

www.tricolbiomedical.com | info@tricolbiomedical.com

Your feedback is important to us.

Please visit www.tricolbiomedical.com/feedback.

