

INTERNAL USE ONLY
— DO NOT PRINT THIS PAGE —

HemCon[®] Dental Dressing **PRO**

10 mm x 12 mm
1 in x 3 in (25 mm x 75 mm)



Manufactured for:

Tricol Biomedical, Inc.

720 SW Washington Street, Suite 200, Portland, OR 97205 U.S.A.

+1.503.245.0459 | Toll Free (U.S. & Canada) 877.247.0196

www.tricolbiomedical.com | info@tricolbiomedical.com

Your feedback is important to us.

Please visit www.tricolbiomedical.com/feedback.



 HemCon
by **Tricol** BIOMEDICAL INC

RM7290 Rev. 1

HemCon[®] Dental Dressing **PRO**

10 mm x 12 mm
1 in x 3 in (25 mm x 75 mm)



Manufactured for:

Tricol Biomedical, Inc.

720 SW Washington Street, Suite 200, Portland, OR 97205 U.S.A.

+1.503.245.0459 | Toll Free (U.S. & Canada) 877.247.0196

www.tricolbiomedical.com | info@tricolbiomedical.com

Your feedback is important to us.

Please visit www.tricolbiomedical.com/feedback.



 HemCon
by **Tricol** BIOMEDICAL INC

RM7290 Rev. 1

HemCon® Dental Dressing PRO

10 mm x 12 mm
1 in x 3 in (25 mm x 75 mm)

Instructions for Use:

1. The HemCon Dental Dressing can be trimmed to fit any bleeding oral wound; dressing should fit loosely in extraction sockets. Both sides of the dressing are active.
2. The dressing is most effective in the presence of enough blood to wet the contact surface of the dressing.
3. Using clean dry instrument or gloved finger, place dressing over the wound site or into the extraction socket. Do not force the dressing into the apex of the socket. The top of the dressing should be flush with the crestal gingiva. (Note: dressing will stick to a wet instrument or glove).
4. Dressing becomes mildly adhesive in contact with blood, apply the dressing once and do not (re-)move, or the dressing will not stick again.
5. Apply pressure to the dressing, either with a clean dry instrument or dry gloved finger, until wound is stabilized (usually 1–2 minutes).
6. Remove pressure with a lateral sweeping motion to avoid dislodging the dressing.
7. Before discharging patient, visually inspect the site to confirm correct positioning of the dressing. For periodontal procedure, a periodontal dressing can be placed over the dressing for protection. Ensure a portion of the dressing is left exposed to oral environment for normal dissolution.
8. The dressing usually dissolves within 48 hours, but may take up to seven days. At seven days, the wound site can be irrigated to ensure removal of any residual material.

Intended Use / Indications: The HemCon Dental Dressing PRO is an oral wound dressing intended as a physical barrier for temporary protection of oral mucosal tissue and to provide pain relief.

Description: The HemCon Dental Dressing PRO is a chitosan-based oral wound dressing intended as a physical barrier for temporary protection of oral mucosal tissue and to provide pain relief.

Cautions / Warnings: Sutures may not be required when using the HemCon Dental Dressing, but can be used at the clinician's discretion. Never suture HemCon Dental Dressing PRO in place using primary closure. An open wound is required when using this device. Do not use if seal is broken and/or contents are wet. Use promptly after opening. Contains chitosan from shellfish. Single use only. Do not resterilize. Not made with natural rubber latex.

Storage Conditions: Best stored at room temperature.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

HemCon® Dental Dressing PRO

10 mm x 12 mm
1 in x 3 in (25 mm x 75 mm)

Instructions for Use:

1. The HemCon Dental Dressing can be trimmed to fit any bleeding oral wound; dressing should fit loosely in extraction sockets. Both sides of the dressing are active.
2. The dressing is most effective in the presence of enough blood to wet the contact surface of the dressing.
3. Using clean dry instrument or gloved finger, place dressing over the wound site or into the extraction socket. Do not force the dressing into the apex of the socket. The top of the dressing should be flush with the crestal gingiva. (Note: dressing will stick to a wet instrument or glove).
4. Dressing becomes mildly adhesive in contact with blood, apply the dressing once and do not (re-)move, or the dressing will not stick again.
5. Apply pressure to the dressing, either with a clean dry instrument or dry gloved finger, until wound is stabilized (usually 1–2 minutes).
6. Remove pressure with a lateral sweeping motion to avoid dislodging the dressing.
7. Before discharging patient, visually inspect the site to confirm correct positioning of the dressing. For periodontal procedure, a periodontal dressing can be placed over the dressing for protection. Ensure a portion of the dressing is left exposed to oral environment for normal dissolution.
8. The dressing usually dissolves within 48 hours, but may take up to seven days. At seven days, the wound site can be irrigated to ensure removal of any residual material.

Intended Use / Indications: The HemCon Dental Dressing PRO is an oral wound dressing intended as a physical barrier for temporary protection of oral mucosal tissue and to provide pain relief.

Description: The HemCon Dental Dressing PRO is a chitosan-based oral wound dressing intended as a physical barrier for temporary protection of oral mucosal tissue and to provide pain relief.

Cautions / Warnings: Sutures may not be required when using the HemCon Dental Dressing, but can be used at the clinician's discretion. Never suture HemCon Dental Dressing PRO in place using primary closure. An open wound is required when using this device. Do not use if seal is broken and/or contents are wet. Use promptly after opening. Contains chitosan from shellfish. Single use only. Do not resterilize. Not made with natural rubber latex.

Storage Conditions: Best stored at room temperature.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.